

A. gossypii Verification PCR

Dirty genomic DNA isolation

- Add 0.2 g of 0.5 mm glass beads into a 1.5 ml tube
- Add 200 μ l DNA extraction buffer
- Add a two pinhead big piece mycelium from the border of a colony and vortex quickly
- Add 100 μ l of water saturated phenol/chloroform
- Vortex for 10 min @ 4 °C
- Centrifuge @ 17000 x g (i.e. 13000 rpm in Eppendorf 5417 R) for 15 min @ RT
- Carefully take 10 μ l of the upper aqueous phase and dilute with 90 μ l dH₂O

DNA Extraction Buffer

- 50 mM NaCl, 1 mM EDTA pH 8, 10 mM Tris · HCl pH 8, 0.5 % Triton X-100
 - 2.922 g NaCl (M_r 58.44 g/mol)
 - 2 ml 0.5 M EDTA pH 8
 - 10 ml 1 M Tris HCl pH 8
 - 50 ml 10 % Triton X-100
- Fill up with dH₂O to 1000 ml
- Autoclave
 - Store @ RT

PCR reaction

- 1 μ l primer 1 (pmol/ μ l)
- 1 μ l primer 2 (pmol/ μ l)
- 5 μ l “dirty” genomic DNA isolate
- 5 μ l dNTPs (2mM each)
- 5 μ l 10x PCR Buffer
- 3 μ l MgCl₂ (25 mM)
- 15 μ l Betaine 5 M
- 0.5 μ l *Taq* DNA polymerase 5 U/ μ l
- 14.5 μ l dH₂O

PCR Cycle

- 95 °C for 5 min
- 95 °C for 45 s
- 51 °C + 0.3 °C per cycle for 2 min
- 72 °C for 1 min
- Repeat steps 2 – 4 34x
- 4 °C unlimited time